

Anti-TBK1 antibody (100-200 aa) (STJ11106133)

STJ11106133

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:2000 ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	56480
Gene Symbol	Tbk1
Uniprot ID	TBK1_MOUSE
Immunogen	
Immunogen Region	100-200 aa
Specificity	A synthetic peptide corresponding to a sequence within amino acids 100-200 of mouse TBK1 (NP_062760.3).
Immunogen Sequence	EPSNAYGLPESEFLIVLRDV VGGMNHLRENGIVHRDIKPG NIMRVIGEDGQSVYKLTDFG AARELEDDEQFVSLYGTEEY LHPDMYERAVLRKDHQKKYG A