

Anti-FAT1 antibody (2700-2900 aa) (STJ11106135)

STJ11106135

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:2000
Range	ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	2195
Gene Symbol	FAT1
Uniprot ID	FAT1_HUMAN
Immunogen	
Immunogen Region	2700-2900 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 2700-2900 of human FAT1 (NP_005236.2).
Immunogen Sequence	LPKFSEPFYFTVSEDVPIG TEIDLIRAEHSGTVLYSLVK GNTPESNRDESFVIDRQSGR LKLEKSLDHETTKWYQFSIL ARCTQDDHEMVASVDVSIQV KDANDNSPVFESSPYEAFIV ENLPGGSRVIQIRASDADSG TNGQVMYSLDQSQSVEVIES FAINMETGWITTLKELDHEK RDNYQIKVVASDHGEKQLS S