

Anti-TRDN antibody (1-101 aa) (STJ11106144)

STJ11106144

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:2000 ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	

TARGET INFORMATION

Gene ID	10345
Gene Symbol	TRDN
Uniprot ID	TRDN_HUMAN
Immunogen	
Immunogen Region	1-101 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1-101 of human TRDN (NP_001242949.1).
Immunogen Sequence	MTEITAEGNASTTTTVIDSK NGSVPKSPGKVLKRTVTEDI VTFSSPAAWLLVIALIITW SAVAIVMFDLVDYKNFSASS IAKIGSDPLKLVRDAMEETT D