

Anti-DDX27 antibody (587-796 aa) (STJ11106148)

STJ11106148

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	ELISA
Host/Source	Rabbit
Reactivity	Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	ELISA: Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	55661
Gene Symbol	DDX27
Uniprot ID	DDX27_HUMAN
Immunogen	
Immunogen Region	587-796 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 587-796 of human DDX27 (NP_060365.7).
Immunogen Sequence	KARILPQDVILKFRDKIEKM EKDVYAVLQLEAEKEMQQS EAQINTAKRLLEKGKEAVVQ EPERSWFQTKEERKKKEKIAK ALQEFDLALRGKKKKRKKFMK DAKKKGEMTAEERSQFEILK AQMFAERLAKRNRRAKRARA MPEEEPVRGPAKKQKQGKKS VFDEELNTSKKALKQYRAG PSFEERKQLGLPHQRRGGNF KSKSRYKRRK