

Anti-TFB2M antibody (20-396 aa) (STJ11106149)

STJ11106149

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:2000
Range	ELISA:Recommended starting concentration is 1 μ g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	64216
Gene Symbol	TFB2M
Uniprot ID	TFB2M_HUMAN
Immunogen	
Immunogen Region	20-396 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 20-396 of human TFB2M (NP_071761.1).
Immunogen Sequence	AGRFCILGSEAAATRKHLPAR NHCGLSDSSPQLWPEPDFRN PPRKASKASLDFKRYVTDRR LAETLAQIYLGKPSRPPHLL LECNPGPGILTQALLEAGAK VVALES DKTFIPHLESLGKN LDGKL RVIHCDFFKLDPRSG GVIKPPAMSSRGLFKNLGIE AVPWTADIP LKVVG MFPSRG EKRALWK LAYDLYSCTSIYK FGRIEVNMFIGEKEFQKLMA DPGNPDLYHVL SVIWLQAC