

Anti-SEMA6D antibody (20-110 aa) (STJ11106150)

STJ11106150

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:2000 ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	

TARGET INFORMATION

Gene ID	80031
Gene Symbol	SEMA6D
Uniprot ID	SEM6D_HUMAN
Immunogen	
Immunogen Region	20-110 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 20-110 of human SEMA6D (NP_705871.1).
Immunogen Sequence	AVSFPEDDEPLNTVDYHYSR QYPVFRGRPSGNSQHRDLF QLMLKIRDTLYIAGRDQVYT VNLNEMPKTEVIPNKLTWR SRQQDRENCAM

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081