

Anti-DDX53 antibody (1-160 aa) (STJ11106152)

STJ11106152

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	IHC-P/ELISA
Host/Source	Rabbit
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	IHC-P:1:50-1:200
Range	ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	168400
Gene Symbol	DDX53
Uniprot ID	DDX53_HUMAN
Immunogen	
Immunogen Region	1-160 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1-160 of human DDX53 (NP_874358.2).
Immunogen Sequence	MSHWAPEWKRAEANPRDLGA SWDVRGSRGSGWSPFGHQG PRAAGSREPPLCFKIKNNMV GVVIGYSGSKIKDLQHSTNT KIQINGESEAKVRIFGNRE MKAKAKAAIETLRKQESYN SESSVDNAASQTPIGRNLGR NDIVGEAEPLSNWDRIRAAV