

Anti-PAM antibody (800 aa) (STJ11106154)

STJ11106154

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:2000 ELISA:Recommended starting concentration is 1 μ g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	

TARGET INFORMATION

Gene ID	5066
Gene Symbol	PAM
Uniprot ID	AMD_HUMAN
Immunogen	
Immunogen Region	800 aa
Specificity	A synthetic peptide corresponding to a sequence within amino acids 800 to the C-terminus of human PAM (NP_620176.1).
Immunogen Sequence	GRVLGRFRGKGSGGLNLGNF FASRKGYSRKGFDRLSTEGS DQEKEDDGSESEEEYSAPLP ALAPSSS