

Anti-IGF2BP1 antibody (478-577 aa) (STJ11106160)

STJ11106160

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/IHC-P/IF/ICC/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:2000 IHC-P:1:50-1:100 IF/ICC:1:50-1:200 ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS, pH 7.3.
Isotype	IgG
Storage Instruction	

TARGET INFORMATION

Gene ID	10642
Gene Symbol	IGF2BP1
Uniprot ID	IF2B1_HUMAN
Immunogen	
Immunogen Region	478-577 aa
Specificity	A synthetic peptide corresponding to a sequence within amino acids 478-577 of human IGF2BP1 (NP_006537.3).
Immunogen Sequence	NFFGPKKEEVKLETHIRVPAS AAGRVIKGGKTVNELQNLT AAEVVVPRDQTPDENDQVIV KIIGHFYASQMAQRKIRDIL AQVKQQHQKQSNQAQARRK