

Anti-AEN antibody (1-325 aa) (STJ11106161)

STJ11106161

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:2000
Range	ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	64782
Gene Symbol	AEN
Uniprot ID	AEN_HUMAN
Immunogen	
Immunogen Region	1-325 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1-325 of human AEN (NP_073604.3).
Immunogen Sequence	MVPREAPESAQCCLPSLTIP NAKDVLRLKRRSRQHQRF MARKALLQEQGLLSMPPEPG SSPLPTPFGAATATEAASSG KQCLRAGSGSAPCSRPPAPG KASGPLPSKCAIDCEMVGT GPRGRVSELARCSIVSYHGN VLYDKYIRPEMPIADYRTRW SGITRQHMRKAVPFQVAQKE ILKLLKGKVVVGHALHNDFQ ALKYVHPRSQTRDITTYVPNF LSEPLHTRARVSLKDLAL