

## Anti-VPS26B antibody (1-160 aa) (STJ11106162)

STJ11106162

### GENERAL INFORMATION

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	
<b>Applications</b>	ELISA
<b>Host/Source</b>	Rabbit
<b>Reactivity</b>	Human

### PRODUCT PROPERTIES

<b>Clonality</b>	Polyclonal
<b>Clone ID</b>	
<b>Concentration</b>	Lot specific
<b>Conjugation</b>	Unconjugated
<b>Purification</b>	Affinity purification
<b>Dilution Range</b>	ELISA: Recommended starting concentration is 1 $\mu$ g/mL. Please optimize the concentration based on your specific assay requirements.
<b>Formulation</b>	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
<b>Isotype</b>	IgG
<b>Storage Instruction</b>	

### TARGET INFORMATION

<b>Gene ID</b>	112936
<b>Gene Symbol</b>	VPS26B
<b>Uniprot ID</b>	VP26B_HUMAN
<b>Immunogen</b>	
<b>Immunogen Region</b>	1-160 aa
<b>Specificity</b>	Recombinant fusion protein containing a sequence corresponding to amino acids 1-160 of human VPS26B (NP_443107.1).
<b>Immunogen Sequence</b>	MSFFFGQGSVEVEILLNDAE SRKRAEHTEDGKKEKYFLF YDGETVSGKVSLALKNPNKR LEHQGIKIEFIGQIELYYDR GNHHEFVSLVKDLARPGEIT QSQAQDFEFTHVEKPYESYT GQNVKRLRYFLRATISRRLND VKEMDIVVHTLSTYPELNS