

Anti-NAP1L4 antibody (1-375 aa) (STJ11106165)

STJ11106165

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB: 1:500-1:2000
Range	ELISA: Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	4676
Gene Symbol	NAP1L4
Uniprot ID	NP1L4_HUMAN
Immunogen	
Immunogen Region	1-375 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1-375 of human NAP1L4 (NP_005960.1).
Immunogen Sequence	MADHSFSDGVPDSVEAAKN ASNTEKLTQVMQNPRVLAA LQERLDNVPHTPSSYIETLP KAVKRRINALKQLQVRCAHI EAKFYEEVHDLERKYAALYQ PLFDKRRREFITGDVEPTDAE SEWHSENEKEELAGDMKSK VVTEKAAATAEEPDPKGIP EFWFTIFRNVDMLSLVQEY DEPIKHLQDIKVFSDPGQ PMSFVLEFHFEPNDYFTNSV LTKTYKMKSEPDKADPFSF