

Anti-LMTK2 antibody (1100-1300 aa) (STJ11106170)

STJ11106170

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:200-1:2000
Range	ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	22853
Gene Symbol	LMTK2
Uniprot ID	LMTK2_HUMAN
Immunogen	
Immunogen Region	1100-1300 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1100-1300 of human LMTK2 (NP_055731.2).
Immunogen Sequence	YRDSAYFSDNDSEPEKRSEE VPGTSPSALVLVQEQPLPEP VLPEQSPAAQDSCLEARKSQ PDESCLSALHNSSDLELRAT PEPAQTGVPQQVHPTDEAS SPWSVLNAELSSGDDFETQD DRPCTLASTGTNTNELLAYT NSALDKSLSSHSEGPKLKEP DIEGKYLGKLGVSGLDLSE DGMDADEEDENSDDSDLELR A