

Anti-SHC3 antibody (310-450 aa) (STJ11106171)

STJ11106171

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:200-1:2000
Range	ELISA:Recommended starting concentration is 1 μ g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	53358
Gene Symbol	SHC3
Uniprot ID	SHC3_HUMAN
Immunogen	
Immunogen Region	310-450 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 310-450 of human SHC3 (NP_058544.3).
Immunogen Sequence	QCPTKIPALHDRMQSLDEPW TEEEGDGS DHPYNSIPSKM PPPGGFLDTRLKPRPHAPDT AQFAGKEQTYQGRHLGDTF GEDWQQTPLRQGSSDIYSTP EGKLHVAPTGEAPTYVNTQQ IPPQAWPAAVSSAESSPRKD L