

## Anti-XPO4 antibody (400-700 aa) (STJ11106174)

STJ11106174

### GENERAL INFORMATION

|                          |                    |
|--------------------------|--------------------|
| <b>Product Type</b>      | Primary antibodies |
| <b>Short Description</b> |                    |
| <b>Applications</b>      | WB/ELISA           |
| <b>Host/Source</b>       | Rabbit             |
| <b>Reactivity</b>        | Human              |

### PRODUCT PROPERTIES

|                      |                                                                                                                                   |
|----------------------|-----------------------------------------------------------------------------------------------------------------------------------|
| <b>Clonality</b>     | Polyclonal                                                                                                                        |
| <b>Clone ID</b>      |                                                                                                                                   |
| <b>Concentration</b> | Lot specific                                                                                                                      |
| <b>Conjugation</b>   | Unconjugated                                                                                                                      |
| <b>Purification</b>  | Affinity purification                                                                                                             |
| <b>Dilution</b>      | WB:1:200-1:2000                                                                                                                   |
| <b>Range</b>         | ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements. |
| <b>Formulation</b>   | PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.                                                                                  |
| <b>Isotype</b>       | IgG                                                                                                                               |
| <b>Storage</b>       |                                                                                                                                   |
| <b>Instruction</b>   |                                                                                                                                   |

### TARGET INFORMATION

|                           |                                                                                                                                                                                                                                                             |
|---------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Gene ID</b>            | 64328                                                                                                                                                                                                                                                       |
| <b>Gene Symbol</b>        | XPO4                                                                                                                                                                                                                                                        |
| <b>Uniprot ID</b>         | XPO4_HUMAN                                                                                                                                                                                                                                                  |
| <b>Immunogen</b>          |                                                                                                                                                                                                                                                             |
| <b>Immunogen Region</b>   | 400-700 aa                                                                                                                                                                                                                                                  |
| <b>Specificity</b>        | Recombinant fusion protein containing a sequence corresponding to amino acids 400-700 of human XPO4 (NP_071904.4).                                                                                                                                          |
| <b>Immunogen Sequence</b> | MEAYDKLLESWLTQVDDKH FHKGFFTQHAVQVFNSYIQC HLAAPDGTRNLTANGVASRE EEEISELQEDDRDQFSDQLA<br>SVGMLGRIAAEHCIPLLTSL LEERVTRLHGQLQRHQQL ASPGSSTVDNKMLDLYEDI HWLILVTGYLLADDTQGETP<br>LIPPEIMEYSIKHSSEVDIN TTLQLGSPGEKASSIPGYN RTDSVIRLLSAILRVSEVES RAIRADLTHLLSPQMGKDI |