

Anti-PARP9 antibody (400-700 aa) (STJ11106176)
STJ11106176

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:200-1:2000
Range	ELISA:Recommended starting concentration is 1 μ g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	83666
Gene Symbol	PARP9
Uniprot ID	PARP9_HUMAN
Immunogen	
Immunogen Region	400-700 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 400-700 of human PARP9 (NP_001139577.1).
Immunogen Sequence	MEIKKETAAEILFDEVLTFADHVKHQLTVKFVIFPTDLEIYKAFSSEMAKRSKMLSLNNYSVPQSTREEKRENGLEARS PAINLMGFNVEEMYEAHAWIQRILSLQNHIIENNHILYLGRKEHDILSQLQKTSSVSIT EIISPGRTELEIEGARADLI EVMNIEDMLCKVQEEMARKKERGLWRSLGQWTIQQQKTQDEMKENIIFLKCPV/PPTQELLDQKKQFEKCGQLQVLKVEK