

Anti-ZNF564 antibody (30-110 aa) (STJ11106178)

STJ11106178

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:200-1:2000 ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	

TARGET INFORMATION

Gene ID	163050
Gene Symbol	ZNF564
Uniprot ID	ZN564_HUMAN
Immunogen	
Immunogen Region	30-110 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 30-110 of human ZNF564 (NP_659413.1).
Immunogen Sequence	RDVMRETFRNLCVGGKKWED QSIEDWYKNQGRILRNHMEELGSESKEYDQCGEAFSQILN LNLNKKIPTIVRPECSCSG K