

Anti-DAGLB antibody (200-400 aa) (STJ11106179)

STJ11106179

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:200-1:2000
Range	ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	221955
Gene Symbol	DAGLB
Uniprot ID	DGLB_HUMAN
Immunogen	
Immunogen Region	200-400 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 200-400 of human DAGLB (NP_631918.3).
Immunogen Sequence	LCCIGKDDHTRVAFSSTAE LFSTYFSDTDLVPSDIAAGL ALLHQQQDNIRNNQEPQV CHAPGSSQEADLDAELENCH HYMQFAAAAYGWPLYIYRNP LTGLCRIGGDCCRSRTTDYD LVGGDQLNCHFGSILHTTGL QYRDFIHVSFHDKVYELPFL VALDHRKESVVAVRGTM SL QDVLTDLSAESEVL DVECEV Q