

Anti-PRR11 antibody (250-350 aa) (STJ11106180)

STJ11106180

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:2000 ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	55771
Gene Symbol	PRR11
Uniprot ID	PRR11_HUMAN
Immunogen	
Immunogen Region	250-350 aa
Specificity	A synthetic peptide corresponding to a sequence within amino acids 250-350 of human PRR11 (NP_060774.2).
Immunogen Sequence	IPSPKARNPLVTVSDLQHVT LKPN SKVLSTRVTNVLITPG KSQMDLRKLLRKVDVERSPG GTPLTNKENMETGTGLTPVM TQALRRKFLAHPRSPTPTL P