

Anti-ZBTB18 antibody (140-360 aa) (STJ11106181)

STJ11106181

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:1000
Range	ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	10472
Gene Symbol	ZBTB18
Uniprot ID	ZBT18_HUMAN
Immunogen	
Immunogen Region	140-360 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 140-360 of human ZBTB18 (NP_006343.2).
Immunogen Sequence	KVESLSDGSSHIAGDLPSE DEGEDEKLNILPSKRDAAE PGNMWMRLPSDSAGIPQAGG EAEPHATAAGKTVASPCSST ESLSQRSVTSVRDSADVDCV LDLSVKSSLSGVENLNSSYF SSQDVLRSNLVQVKVEKEAS CDESVDGTNDYDMEHSTVKE SVSTNNRVQYEPAPLAPLRE DSVLRELDREDKASDDEMMT PESERVQVEGGMESLLPYV S