

Anti-ACAP1 antibody (20-230 aa) (STJ11106182)

STJ11106182

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:2000
Range	ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	9744
Gene Symbol	ACAP1
Uniprot ID	ACAP1_HUMAN
Immunogen	
Immunogen Region	20-230 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 20-230 of human ACAP1 (NP_055531.1).
Immunogen Sequence	SIELVEAEVSELETRLEKLL KLGTGLLESGRHYLAASRAF VVGICDLARLGPPEPMMMAEC LEKFTVSLNHHKLDSHAELLD ATQHTLQQIQTLVKEGLRG FREARRDFWRGAESLEAALT HNAEVPRRRAQEAEAGAAL RTARAGYRGRALDYALQINV IEDKRFDIIMEFVLRLEAQ ATHFQQGHEELSRLSQYRKE LGAQLHQLVLN