

Anti-TRAF3IP3 antibody (290-530 aa) (STJ11106184)
STJ11106184

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:2000
Range	ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	80342
Gene Symbol	TRAF3IP3
Uniprot ID	T3JAM_HUMAN
Immunogen	
Immunogen Region	290-530 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 290-530 of human TRAF3IP3 (NP_079504.2).
Immunogen Sequence	VLEEKMNAEQQLQSTQRLA LAEQKCEEWRSQYEALKEDW RTLGTQHRELSQLHVLQSK LQGADSRDLQMNQALRFLN EHQQLQAKIECLQGDRDLCS LDTQDLQDQLKRSEAEKLT VTRVQQLQGLLQNGSLQLQE QEKLLTKKDQALPVWSPKSF PNEVEPEGTGKEKDWDLRDQ LQKKTLLQLQAKEKECRELHS ELDNLSDEYLSCLRKLQHCR EELNQSQQLPPRRQCGRWL