

Anti-PHF2 antibody (850-950 aa) (STJ11106191)

STJ11106191

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	ELISA
Host/Source	Rabbit
Reactivity	Human/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	ELISA: Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	5253
Gene Symbol	PHF2
Uniprot ID	PHF2_HUMAN
Immunogen	
Immunogen Region	850-950 aa
Specificity	A synthetic peptide corresponding to a sequence within amino acids 850-950 of human PHF2 (NP_005383.3).
Immunogen Sequence	AKNSVDLDDYEEQDHLDAC FKDSYVYPSLESDENPIF KRSRKKRKGSDDPYSPTAR VGPSVPRQDRPVREGTRVAS IETGLAAAAAKLSQQEEQKS K