

Anti-XPR1 antibody (464-631 aa) (STJ11106194)

STJ11106194

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:2000
Range	ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	9213
Gene Symbol	XPR1
Uniprot ID	S53A1_HUMAN
Immunogen	
Immunogen Region	464-631 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 464-631 of human XPR1 (NP_001129141.1).
Immunogen Sequence	DLKMDWGLFDKNAGENTFLR EEIVYPQKAYYYCAIIEDVI LRFQWTQISITSTLLPHS GDIIATVFAPLEVFRRFVWN FFRLENEHLNNGEFAVRD ISVAPLNADDQTLLEQMMDQ DDGVRNRQKNRSWKYNQSIG LRRPRLASQSKARDTKVLIE DTDDEANT