

Anti-BES1 antibody (89 aa) (STJ11106195)

STJ11106195

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	A.thaliana/Rice

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:2000
Range	ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	
Immunogen	89 aa
Region	
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 89-335 of Arabidopsis thaliana BES1 (NP_001077562.1).
Immunogen Sequence	GHKPLPGDMAGSSSRATPYS SHNQSPLSSTFDSPILSYQV SPSSSFSPSPRVGDPHNIS TIFPFLRNGGIPSSLPLRI SNSAPVTTPPVSSPTSRNPKP LPTWESFTKQSMMAAKQSM TSLNYPFYAVSAPASPTHHR QFHAPATIECEDSDSSTVD SGHWISFQKFAQQPFASAM VPTSPTFNLVKPAPQQLSPN TAAIQEIGQSSEFKFENSQV KPWEGERIHDVAMEDLELT