

Anti-PIK3R5 antibody (500-600 aa) (STJ11106196)

STJ11106196

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:1000 ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	23533
Gene Symbol	PIK3R5
Uniprot ID	PI3R5_HUMAN
Immunogen	
Immunogen Region	500-600 aa
Specificity	A synthetic peptide corresponding to a sequence within amino acids 500-600 of human PIK3R5 (NP_055123.2).
Immunogen Sequence	QRRRPFLSGDEDPKASTLRV VVFGSDRISGKVARAYSNLR RLENNRPLLTRFFKLQFFYV PVKRSHGTSPGACPPPRSQT PSPPTDSPRHASPGELGTTP W