

Anti-SV2C antibody (460-580 aa) (STJ11106197)

STJ11106197

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:2000 ELISA:Recommended starting concentration is 1 μ g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	

TARGET INFORMATION

Gene ID	22987
Gene Symbol	SV2C
Uniprot ID	SV2C_HUMAN
Immunogen	
Immunogen Region	460-580 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 460-580 of human SV2C (NP_055794.3).
Immunogen Sequence	PLQSDHEYALLTRNVERDKYA NFTINFTMENQIHTGMEYDN GRFIGVKFKSVTFKDSVFKS CTFEDVTSVNTYFKNCTFID TVFDNTDFEPYKFIDSEFKN CSFFHNKTGCQITFDDDYSA Y