

Anti-Cytokeratin 1 antibody (200-450 aa) (STJ11106198)

STJ11106198

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	

TARGET INFORMATION

Gene ID	3848
Gene Symbol	KRT1
Uniprot ID	K2C1_HUMAN
Immunogen	
Immunogen Region	200-450 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 200-450 of human Cytokeratin 1 (Cytokeratin 1 (KRT1)) (NP_006112.3).
Immunogen Sequence	FLEQQNQVLQTKWELLQQVD TSTRTHNLEPYFESFINNLR RRVDQLKSDQSRDSELKNM QDMVEDYRNKYEDEINKRTN AENEFVTIKKDVGDGAYMTKV DLQAKLDNLQQEIDFLTALY QAELSQMQTQISETNVILSM DNNRSLDLSIIAEVKAQYE DIAQSKAEAESLYQSKYEE LQITAGRHGDSVRNSKIEIS ELNRVIQRLRSEIDNVKKQI SNLQQSISDAEQRGENALK