

Anti-MARCH7 antibody (425-704 aa) (STJ11106200)
STJ11106200

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:2000
Range	ELISA:Recommended starting concentration is 1 μ g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	64844
Gene Symbol	MARCH7
Uniprot ID	MARH7_HUMAN
Immunogen	
Immunogen Region	425-704 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 425-704 of human MARCH7 (NP_073737.1).
Immunogen Sequence	RRRESNEVVHLEAQNDPLGAA ANRPQASAAASSSATTGGSTS DSAQGGGRNTGISGILPGSLF RFAVPPALGSNLTDNVMITV DIIPSGWNSADGKSDTKSA PSRDPERLQKIKESLLLEDS EEEEGDLCRICQMAAASSN LLIEPCKCTGSLQYVHQDCM KKWLQAKINSGSSLEAVTTC ELCKEKLELNLEDFDIHELH RAHANEQAEYEFISSGLYLV VLLHLCEQSFSDMMGNTNE