

Anti-ADRA1D antibody (450-550 aa) (STJ11106202)

STJ11106202

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:2000 ELISA:Recommended starting concentration is 1 μ g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	146
Gene Symbol	ADRA1D
Uniprot ID	ADA1D_HUMAN
Immunogen	
Immunogen Region	450-550 aa
Specificity	A synthetic peptide corresponding to a sequence within amino acids 450-550 of human ADRA1D (NP_000669.1).
Immunogen Sequence	APSSGDAPPGAPLALTALPD PDPEPPGTPEMQAPVASRRK PPSAFREWRLGPFRRPTTQ LRAKVSSLSHKIRAGGAQRA EAACAQRSEVEAVSLGVPHE V