

Anti-COL4A5 antibody (1462-1691 aa) (STJ11106204)

STJ11106204

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	ELISA
Host/Source	Rabbit
Reactivity	Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	ELISA: Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	

TARGET INFORMATION

Gene ID	1287
Gene Symbol	COL4A5
Uniprot ID	CO4A5_HUMAN
Immunogen Region	1462-1691 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1462-1691 of human COL4A5 (NP_203699.1).
Immunogen Sequence	SSVAHGFLITRHSQTDDAPQ CPQGTQVYEGFSLLYVQGN KRAHGQDLGTAGSCLRRFST MPFMFCNINNVCNFASTRNDY SYWLSTPEPMPMSMQPLKGQ SIQPFISRCAVCEAPAVIA VHSQTIQIPHCPQGWDSLWI GYSFMMHTSAGAECSGQALA SPGSCLEEFRSAPFIECHGR GTCNYYANSYSFWLATVDVS DMFSKPQSETLKAGDLRTRI SRCQVCMKRT