

Anti-DGKB antibody (544-773 aa) (STJ11106205)

STJ11106205

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	IF/ICC/ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	IF/ICC:1:50-1:200
Range	ELISA:Recommended starting concentration is 1 μ g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	1607
Gene Symbol	DGKB
Uniprot ID	DGKB_HUMAN
Immunogen	
Immunogen Region	544-773 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 544-773 of human DGKB (NP_663733.1).
Immunogen Sequence	LMKILKDIENSTEIMLDRWK FEVIPNDKDEKGDVPYSII NNYFSIGVDASIAHRFHIMR EKHPEKFNSRMKNKFWYFEF GTSETFSATCKKLHESVEIE CDGVQIDLINISLEGIAILN IPSMHGGSNLWGESKKRRSH RRIEKKGSDKRTTVDAKEL KFASQDLSDQLLEVVGLEGA MEMGQIYTGLKSAGRRLAQC SCVVIRTSKSLPMQIDGEPW MQTPCTVSTE