

Anti-FOXA3 antibody (1-100 aa) (STJ11106206)

STJ11106206

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:2000 ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	3171
Gene Symbol	FOXA3
Uniprot ID	FOXA3_HUMAN
Immunogen	
Immunogen Region	1-100 aa
Specificity	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human FOXA3 (NP_004488.2).
Immunogen Sequence	MLGSVKMEAHDLAEWSYYPE AGEVYSPVTPVPTMAPLNSY MTLNPLSSPYPPGGLPASPL PSGPLAPPAPAAPLGPTFPGLGVSGGSSSSGYGAPGPGLV