

Anti-ZFYVE9 antibody (1-260 aa) (STJ11106209)

STJ11106209

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	ELISA: Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	9372
Gene Symbol	ZFYVE9
Uniprot ID	ZFYV9_HUMAN
Immunogen	
Immunogen Region	1-260 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1-260 of human ZFYVE9 (NP_004790.2).
Immunogen Sequence	MENYFQAEAYNLDKVLDEFE QNEDETVSSTLLDTKWNKIL DPPSHRLSFNPTLASVNESA VSNESQPQLKVFLAHSAPL TTEEDHCANGQDCNLNPEI ATMWIDENAVAEDQLIKRNY SWDDQCSAVEVGEKKCGNLA CLPDEKNVLVVAVMHNC DKR TLQNDLQDCNNYNSQLMDA FSCSLDNENRQTDQFSFSIN ESTEKMNSEKQMDPLNRPK TEGRSVNHLCPTSSDSLAS