

Anti-HOMER2 antibody (144-343 aa) (STJ11106210)

STJ11106210

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	ELISA: Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	

TARGET INFORMATION

Gene ID	9455
Gene Symbol	HOMER2
Uniprot ID	HOME2_HUMAN
Immunogen Region	144-343 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 144-343 of human HOMER2 (NP_004830.2).
Immunogen Sequence	AGPANTHLKSENDKLIKALT QSAANVKKWEIELQTLRESN ARLTTALQESAASVEQWKRQ FSICRDENDRLRNKIDELEE QCSEINREKEKNTQLKRRIE ELEAELREKETELKDLRKQS EIIPQLMSECEYVSEKLEAA ERDNQNLEDKVRSLKTDIEE SKYRQRHLKVELKSFLEVLGDKIDDLHDFRRGLSKLGTDN