

Anti-Paxillin antibody (54-253 aa) (STJ11106218)

STJ11106218

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:2000
Range	ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	5829
Gene Symbol	PXN
Uniprot ID	PAXI_HUMAN
Immunogen	
Immunogen Region	54-253 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 54-253 of Paxillin (NP_001074324.1).
Immunogen Sequence	SSEALNGTILDPLDQWQPSS SRFIHQQPQSSSPVYGSSAK TSSVSNPQDSVGSPCSRVE EEHVYSFPNKOKSAEPSPTV MSTSLGSLNSELDRLLLELN AVQHNPPGFPADEANSSPPL PGALSPLYGVPETNSPLGGK AGPLTKEKPKRNGGRGLEDV RPSVESLLDELESSVPSPVP AITVNQGEEMSSPQRTSTQQ