

Anti-KISS1 antibody (1-138 aa) (STJ11106223)

STJ11106223

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	ELISA
Host/Source	Rabbit
Reactivity	Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	

TARGET INFORMATION

Gene ID	3814
Gene Symbol	KISS1
Uniprot ID	KISS1_HUMAN
Immunogen	
Immunogen Region	1-138 aa
Specificity	A synthetic peptide corresponding to a sequence within amino acids 1-138 of human KISS1 (NP_002247.3).
Immunogen Sequence	MNSLVSWQLLLFLCATHFGE PLEKVASVGNSRPTGQQLES LGLLAPGEQSLPCTERKPAA TARLSRRGTSLSPPPESGSGS
Sequence	PQQPGLSAPHSRQIPAPQGA VLVQREKDLPNYNWNSFGLR FGKREAAPGNHGRSAGRG