

Anti-ARF4 antibody (1-180 aa) (STJ11106231)

STJ11106231

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	ELISA: Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	378
Gene Symbol	ARF4
Uniprot ID	ARF4_HUMAN
Immunogen	
Immunogen Region	1-180 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1-180 of human ARF4 (NP_001651.1).
Immunogen Sequence	MGLTISLFSRLFGKKQMRI LMVGLDAAGKTTILYKLKLG EIVTTIPTIGFNVETVEYKN ICFTVWDVGGQDRIRPLWKH YFQNTQGLIFVDSNDRERI QEVADELQKMLLVDELDAV LLLFANKQDLPNAMAISEMT DKLGLQSLRNRTWYVQATCA TQGTGLYEGLDWLSNELSKR