

Anti-ARF4 antibody (1-180 aa) (STJ11106232)

STJ11106232

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	IF/ICC/ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	IF/ICC:1:50-1:200
Range	ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	378
Gene Symbol	ARF4
Uniprot ID	ARF4_HUMAN
Immunogen	
Immunogen Region	1-180 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1-180 of human ARF4 (NP_001651.1).
Immunogen Sequence	MGLTISSLFSRLFGKKQMRI LMVGLDAAGKTTILYKLKLG EIVTTIPTIGFNVETVEYKN ICFTVWDVGGQDRIRPLWKH YFQNTQGLIFWDSNDRERI QEVADELQKMLLVDELDAV LLLFANKQDLPNAMAISEMT DKLGLQSLRNRTWYVQATCA TQGTGLYEGLDWLSNELSKR