

Anti-C4BPA antibody (298-597 aa) (STJ11106233)

STJ11106233

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/IF/ICC/ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB: 1:500-1:2000
Range	IF/CC: 1:50-1:200
	ELISA: Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	722
Gene Symbol	C4BPA
Uniprot ID	C4BPA_HUMAN
Immunogen	
Immunogen Region	298-597 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 298-597 of human C4BPA (NP_000706.1).
Immunogen Sequence	SCINLPDIPHASWETYPPT KEDVYVVGTVLRYRCHPGYK PTTDEPTTVICQKNLRWTPY QGCEALCCPEPKLNNGEITQ HRKSRPANHCYFYGDEISF SCHETSRFSAICQGDGTWSP RTPSCGDICNFPPKIAHGHI KQSSSYFFKEEIIIECDKG YILVGQAKLSCSYSHWSAPA PQCKALCRKPELVNGRLSVD KDQYVEPENVTIQCDSGYGV VGPQSITCSGNRTWYPEVP