

Anti-CAPN3 antibody (70-309 aa) (STJ11106234)

STJ11106234

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:2000
Range	ELISA:Recommended starting concentration is 1 μ g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	825
Gene Symbol	CAPN3
Uniprot ID	CAN3_HUMAN
Immunogen	
Immunogen Region	70-309 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 70-309 of human CAPN3 (NP_775111.1).
Immunogen Sequence	EEVENTISVDRPVKKKTKP IIFVSDRANSNKELGVDQES EEGKGKTSPOKQKSPQPP GSDDQEEQQFRNIFKQI AGDDMEICADELKKVLNTVV NKHKDLKTHGFTLESCRSMI ALMDTDGSGKLNLFHHLW NKIKAWQKIFKHYDQSGT INSYEMRNAVNDAGFHLNNQ LYDIITMRYADKHMNIDFDS FICCFVRLEGMFRAFHAFDK DGDGIIKLNVLWLQLTMY