

Anti-CCNT2 antibody (250-550 aa) (STJ11106239)

STJ11106239

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:2000
Range	ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	905
Gene Symbol	CCNT2
Uniprot ID	CCNT2_HUMAN
Immunogen	
Immunogen Region	250-550 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 250-550 of human CCNT2 (NP_001232.1).
Immunogen Sequence	RLKKIRNWRANQAARKPKVD GQVSETPLLGSSLVQNSILV DSVTGVPTNPSFQKPSTSAF PAPVPLNSGNISVQDSHTSD NLSMLATGMPSTSYGLSSHQ EWPQHQDSARTEQLYSQKQE TSLSGSQYNINFQQGPSISL HSGLHHRPDKISDHSSVKQE YTHKAGSSKHHGPISTTPGI IPQKMSLDKYREKRKLETLD LDVRDHYIAAQVEQQHKQGQ SQAASSSVTSPIKMKIPI