

Anti-COL15A1 antibody (200-450 aa) (STJ11106243)
STJ11106243

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	IHC-P/ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	IHC-P:1:50-1:200
Range	ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	1306
Gene Symbol	COL15A1
Uniprot ID	COFA1_HUMAN
Immunogen	
Immunogen Region	200-450 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 200-450 of human COL15A1 (NP_001846.3).
Immunogen Sequence	AGIFMGNAGATGLERFTGSL QQLTVHPDPRTPEELCDPEE SSASGETSGLQADGVAEIL EAVTYTQASPKEAKVEPINT PPTPSSPFDEMLSGEPVPE GTLETTNMSIIQHSSPKQGS GEILNDTLEGVHSVVDGDPIT DSGSGAGAFDIAEEKNLAA TAAGLAEVPISTAGEAEASS VPTGGPTLSMSTENPEEGVT PGPDNEERLAATAAGEAEAL ASMPGEVEASGVAPGELDL