

Anti-DGKG antibody (1-250 aa) (STJ11106245)

STJ11106245

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:2000
Range	ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	1608
Gene Symbol	DGKG
Uniprot ID	DGKG_HUMAN
Immunogen	
Immunogen Region	1-250 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1-250 of human DGKG (NP_001074214.1).
Immunogen Sequence	MGEERWVSLTPEEFDQLQKY SEYSSKKIKDALTEFNEGGS LKQYDPHEPISYDVFKLFMR AYLEVDLPOPLSTHLFLAFS QKPRHETSDHPTEGASNSEA NSADTNIQNADNATKADEAC APDTESNMAEKQAPAEQVA ATPLEPPVPRSSSSSPVVY LKDVCYLSLLETGRPQDKL EFMFRLYDSDENGLLDQAEM DCIVNQMLHIAQYLEWDPTL LRPILKEMQLQGMDYDRDGF