

Anti-ENO1 antibody (1-434 aa) (STJ11106252)

STJ11106252

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB: 1:500-1:2000
Range	ELISA: Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	2023
Gene Symbol	ENO1
Uniprot ID	ENOA_HUMAN
Immunogen	
Immunogen Region	1-434 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1-434 of human ENO1 (NP_001419.1).
Immunogen Sequence	MSILKIHAREIFDSRGNPTV EVDLFTSKGLFRAAVPSGAS TGIYEALRLDNDKTRYMGK GVS KAVEHINKTIAPALVSK KLNVTQEKEIDKLMIEMDGT ENKSKFGANAILGVSLAVCK AGAVEKGVPLYRHIADLAGN SEVILPVPAFNVINGGSHAG NKLAMQEFMILPVGAANFRE AMRIGAEVYHNLKNVIKEY GKDATNVGDEGGFAPNILEN KEGLELLKTAIGKAGYTDK