

Anti-GRM6 antibody (480-580 aa) (STJ11106256)

STJ11106256

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	IF/ICC/ELISA
Host/Source	Rabbit
Reactivity	Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	IF/ICC:1:50-1:200 ELISA:Recommended starting concentration is 1 μ g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	2916
Gene Symbol	GRM6
Uniprot ID	GRM6_HUMAN
Immunogen	
Immunogen Region	480-580 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 480-580 of human GRM6 (NP_000834.2).
Immunogen Sequence	GSASSGGYQAVGQWAETLRL DVEALQWSGDPHEVPSSLCS LPCGPGERKKMKVGVPCQWH CEACDGYRFQVDEFTCEACP GDMRPTPNHTGCRPTPVRL S