

Anti-IFNAR2 antibody (27-243 aa) (STJ11106259)

STJ11106259

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	ELISA: Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	

TARGET INFORMATION

Gene ID	3455
Gene Symbol	IFNAR2
Uniprot ID	INAR2_HUMAN
Immunogen Region	27-243 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 27-243 of human IFNAR2 (NP_997467.1).
Immunogen Sequence	ISYDSPDYTDESCTFKISLR NFRSILSWELKNHSIVPTHY TLLYTIMSKPEDLKVVKNCA NTTRSFCDLTDEWRSTHEAY VTVLEGFSGNITLFCSHNF WLAIMSFEPPEFEIVGFTN HINVMVKFPSIVEEELQFDL SLVIEEQSEGIVKKHKPEIK GNMSGNFYIIDKLIPNTNY CVSVYLEHSDEQAVIKSPLK CTLLPPGQESESAESAK