

Anti-IFNAR2 antibody (27-243 aa) (STJ11106260)

STJ11106260

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:2000
Range	ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	3455
Gene Symbol	IFNAR2
Uniprot ID	INAR2_HUMAN
Immunogen	
Immunogen Region	27-243 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 27-243 of human IFNAR2 (NP_997468.1).
Immunogen Sequence	ISYDSPDYTDESCTFKISLR NFRSILSWELKNHSIVPTHY TLLYTIMSKPEDLKVKNCANITRSFCDLTDEWRSTHEAY VTVLEGFSGNTTLFSCSHNF WLAIMSFEPPEFEIVGFTN HINVMVKFPSIVEEELQFDL SLVIEEQSEGIVKKHKPEIK GNMSGNFYIIDKLIPNTNY CVSVYLEHSDEQAVIKSPLK CTLLPPGQESESAESAK