

Anti-ADCY9 antibody (550-780 aa) (STJ11106264)

STJ11106264

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:2000
Range	ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	115
Gene Symbol	ADCY9
Uniprot ID	ADCY9_HUMAN
Immunogen	
Immunogen Region	550-780 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 550-780 of human ADCY9 (NP_001107.2).
Immunogen Sequence	VIERLGQSVVADQLKGLKTY LISGQRAKESRCSCAEALLS GFEVIDGSQVSSGPRGQGTAA SSGNVSDLAQTVKTFDNLKT CPSCGITFAPKSEAGAEGGA PQNGCQDEHKNSTKASGGPN PKTQNGLLSPPQEEKLTNSQ TSLCEILQEKGWAGVSLDQ SALLPLRFKNIREKTDAHFV DVIKEDSLMKDYFFKPPINQ FSLNFLDQELERSYRTSYQE EVIKNSPVKTF