

Anti-PEPD antibody (1-260 aa) (STJ11106268)

STJ11106268

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB: 1:500-1:2000
Range	ELISA: Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	5184
Gene Symbol	PEPD
Uniprot ID	PEPD_HUMAN
Immunogen	
Immunogen Region	1-260 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1-260 of human PEPD (NP_000276.2).
Immunogen Sequence	MAAATGPSFWLGNETLKVPL ALFALNRQRLCERLRKNPAV QAGSIVVLQGGREETQRYCTD TGVLFQESFFHWAFGVTEP GCYGVIDVDTGKSTLFVPRP PASHATWMGKIHSKEHFKEK YAVDDVQYVDEIASVLTQK PSVLLTLRGVNTDSGSVCRE ASFDGISKFEVNNILHPEI VECRVFKTDMELEVLRYTNK ISSEAHREVMKAVKVGMMKEY ELESLEHYCYSRGGMRHS